

HaemoCare Install Guide

HaemoCare Team

May 2026

CONTENTS

HaemoCare Install Guide	
Before the patient arrives	
The install QR code	
Things to tell the patient up front	
Step 1 — Check the phone	
Step 2 — Scan the QR code	
Step 3 — Download the APK	
Step 4 — Allow installation from this source	
Step 5 — Install	
Step 6 — First launch (do this together)	
6.1 — Switch language if needed	
6.2 — Create the account	
6.3 — Complete the medical profile	
6.4 — PDPA consent	
Step 7 — Turn off “install from unknown sources”	
Step 8 — Sanity check together	
Common issues and how to handle them	
When you are done	
Tear-off quick checklist	

HaemoCare Install Guide

For team members helping patients install the app on Android.

This guide walks you, the staff member, through helping a patient install HaemoCare on their Android phone for the first time. Sit beside the patient, follow it end to end. The whole flow takes about five to ten minutes.

WHO DOES WHAT

You handle the install link and any tricky Android dialogs. The patient types their own email and password, fills in their own medical profile, and gives their own PDPA consent. Never type the patient's password for them.

Before the patient arrives

- 1 **Print the install-QR handout** (shown below) or have it open on your own phone. The QR points to the latest HaemoCare release on GitHub, where the APK is attached. The same QR keeps working when new versions ship, so a sheet printed today is still valid next month.
- 2 **Have the patient's email and phone number on hand** (in case you need to reach them).
- 3 **Find a quiet spot with reliable Wi-Fi.** The APK is roughly 60 to 80 MB.
- 4 **Block five to ten minutes** of uninterrupted time.

TIP

Ask the patient to charge their phone above 30% before the session, and to connect to your facility's Wi-Fi when they arrive.

THE INSTALL QR CODE



Scan with the patient's Android camera

github.com/ggmrcodes/TNDH_Project/releases/latest

SOURCE FILE

The master QR lives at `HaemoCare/assets/handout/qr-install.png` (600x600 PNG, B/W). To regenerate or point it at a different mirror (e.g. a hospital intranet), run `bash scripts/regenerate-qr.sh` with the `REGENERATE_QR_URL` environment variable set.

Things to tell the patient up front

- **Android only at this stage.** If the patient has an iPhone, stop here and flag your team lead. There is no iOS install path yet.
 - **The app comes as an APK, not from Google Play.** This is normal during the pilot. We turn the “install from unknown sources” permission ON for one install and OFF again at the end.
 - **The patient creates their own account.** You should not type the patient’s password. They choose it and remember it.
 - **The patient owns their data.** They can export everything as PDF, or delete the account entirely, from inside the app at any time.
-

Step 1 — Check the phone

1 Check the phone meets requirements

Confirm all of the following on the patient's phone:

- Android **8.0 (Oreo) or newer**. Check at **Settings** → **About phone** → **Android version**.
- At least **200 MB free storage**.
- Connected to **Wi-Fi**.

STOP IF

Android version is below 8.0. Flag your team lead, do not proceed.

Step 2 — Scan the QR code

2 Scan the install QR with the patient's camera

Hold up the printed QR sheet (or your phone showing it). Have the patient:

1. Open the **Camera** app on their Android phone.
2. Point it at the QR code, hold steady for a second.
3. Tap the link popup that appears at the bottom of the camera screen.

The link opens **GitHub** in Chrome (or Samsung Internet) and lands on the **latest HaemoCare release** page. The page lists the release notes at the top and the downloadable files lower down.

NO CAMERA POPUP?

Some older camera apps don't auto-detect QR codes. Open the **Google Lens** shortcut from the camera, or download **Google Lens** from Play Store. Lens reads any QR and offers the link.

IF THE PAGE DOESN'T LOAD

Check Wi-Fi. Open `google.com` first to confirm internet is working, then try scanning again.

Step 3 — Download the APK

3 Tap the APK on the release page

On the GitHub release page, scroll down to the **Assets** section. The patient taps the file ending in `.apk` (it will be named something like `haemocare-v1.0.0.apk`).

Chrome starts downloading the file. A small progress notification appears at the top of the screen.

Chrome will likely show:

"This type of file can harm your device. Do you want to keep haemocare.apk anyway?"

This is a generic warning for any APK downloaded outside Play Store. The patient taps **Keep** (or **Download anyway**).

Wait 30 to 90 seconds for the download to finish.

Step 4 — Allow installation from this source

4 Turn ON "install from unknown sources" (one time)

When the download finishes, the patient taps the **haemocare.apk** notification, or opens **Files** → **Downloads** and taps the file.

Android will say something like:

"For your security, your phone currently isn't allowed to install unknown apps from this source."

This is expected. From that dialog:

1. Tap **Settings**.
2. The phone opens **Install unknown apps** for Chrome (or whichever browser was used).
3. Flip the toggle to **Allow from this source** → **ON**.
4. Tap **back** to return to the install prompt.

NOTE FOR STAFF

We turn this permission OFF again in **Step 7**. Leaving it on isn't dangerous, but turning it off after the install is good hygiene.

Step 5 — Install

5 Tap Install, dismiss Play Protect if it appears

The patient taps **Install** on the Android install dialog.

Some phones show a **Play Protect** screen first: *"Play Protect didn't recognize this app."* The patient taps **Install anyway**. This is normal for apps not yet on Play Store; it's a one-time scan.

When installation finishes, tap **Open** to launch HaemoCare. If the patient accidentally taps **Done**, they can launch HaemoCare from the app drawer like any other app.

Step 6 — First launch (do this together)

6 Walk the patient through first launch

The first time HaemoCare opens, the patient sees a splash screen for about a second and a half, then the **Sign In** screen. Help the patient through these four screens in order. The patient does the tapping; you stand beside them.

TH EN

H HaemoCare
ผู้ช่วยการรับเลือดของคุณ

เข้าสู่ระบบ

อีเมล

demo@haemocare.app

รหัสผ่าน

.....

→] เข้าสู่ระบบ

ยังไม่มีบัญชี? สร้างบัญชี

เป็นไปตาม PDPA

6.1 — SWITCH LANGUAGE IF NEEDED

Top right of the sign-in screen has an **EN/TH** toggle. The app defaults to **Thai (TH)**; tap **EN** if the patient prefers English.

The screenshot here shows the Thai version. The app name, the tagline (“ผู้ช่วยการรับเลือดของคุณ”), the email + password fields, and the “ยังไม่มีบัญชี? สร้างบัญชี” link at the bottom (English: “Don’t have an account? Create Account”) are all visible.

6.2 — CREATE THE ACCOUNT

Tap “สร้างบัญชี” at the bottom of the sign-in screen. The signup screen has three fields:

- **อีเมล (Email)** — the patient’s personal email. They will need to remember this.
- **รหัสผ่าน (Password)** — at least 8 characters. *The patient chooses. You do not.*
- **ยืนยันรหัสผ่าน (Confirm Password)** — same value again.

Tap “สร้างบัญชี” (Create Account). The app confirms the account and returns to the sign-in screen. The patient now signs in with the same email and password.

PRIVACY

If the patient is shy about typing in front of you, step back and let them type privately. You do not need to see their password. If they get an “invalid email or password” error after signing up, the most common cause is a trailing space in the email — have them retype carefully.

6.3 — COMPLETE THE MEDICAL PROFILE

After signing in, HaemoCare shows the **Complete Your Profile** screen (Thai: “เสร็จสิ้นโปรไฟล์ของคุณ”). The patient enters:

- **Full Name** — legal name as it appears at the hospital.
- **Blood Type and Rh Factor** — pick from the list.
- **Antibodies on file** — tap **Add antibody** for each. If unknown, leave empty; they can edit later from the **Passport** tab.
- **Known Transfusion Reactions** — short description if any.

- **Current Medications** — list them.
- **Days between transfusions** — defaults to **28**. Adjust to the patient's actual schedule. This is what HaemoCare uses to detect overdue visits and bump symptom severity.

Help the patient here. Most of this is on their hospital paperwork. **If they don't know a field, leave it blank rather than guessing.** They can edit later from **Passport** → **Edit Profile**.

Tap **Complete Profile**.

6.4 — PDPA CONSENT

READ THIS SCREEN WITH THE PATIENT — DO NOT SKIP

The PDPA Data Privacy Consent screen explains what HaemoCare collects, why, who can access it, the patient's PDPA rights, and how long data is retained. The patient must tick the consent checkbox and tap "ฉันยอมรับและดำเนินการต่อ" (I Agree & Continue) before they can use the app.

If the patient does not consent, **do not pressure them**. They can leave the app installed and come back another day. Tell them the app cannot proceed without consent, and that they're also free to choose not to install at all.

Step 7 — Turn off “install from unknown sources”

7 Restore the patient's security setting

Once HaemoCare is installed and running, undo the permission from Step 4.

Path varies by phone brand and Android version, but typically:

Settings → **Apps** → **Special access** → **Install unknown apps** → **Chrome** → **Allow from this source** → **OFF**

This means future random APKs cannot install themselves without the patient deliberately enabling the permission again.

Step 8 — Sanity check together

8 Quick tour of the four tabs

Before the patient leaves, walk through each tab together. This is the part that builds confidence.

Tab	Show the patient
Passport (บัตรประจำตัว)	Blood type, antibodies, Rh, reactions. The Show QR Code button. The Export PDF button.
Log (บันทึก)	The Log Symptoms button. Explain the 72-hour active monitoring window after each transfusion.
Appointments (นัดหมาย)	The Add Appointment button and the Import button (for hospital systems).
History (ประวัติ)	How transfusion records appear here. The Scan option for photographing transfusion slips.

Then find **Privacy Settings** from the Passport tab. Show the patient where **Export My Data** and **Delete Account & Data** live. These are their PDPA rights, one tap each.

Finally, ask the patient to sign in once on their own while you watch, so you know they can repeat it at home.

Common issues and how to handle them

"THE CAMERA WON'T SCAN THE QR CODE"

Some older camera apps don't auto-detect QR. Install **Google Lens** from Play Store, point Lens at the QR. As a fallback, type the URL manually into the browser: `github.com/ggm-rcodes/TNDH_Project/releases/latest`.

"THE GITHUB PAGE WON'T LOAD"

Check Wi-Fi. Open `google.com` first to confirm internet is working. If Wi-Fi is the issue, ask whether the patient is comfortable using mobile data for the install (60 to 80 MB). Otherwise reschedule.

"MY PHONE SAYS 'APP NOT INSTALLED'"

Usually a partial download. Open the **Files** app, delete `haemocare.apk` from Downloads, then redo Step 3.

"I CAN'T FIND THE SETTINGS OPTION AFTER TAPPING THE APK"

Some phones hide it behind a three-dot menu in the dialog. Or open **Settings** → **Apps** → **Special access** → **Install unknown apps** manually, find Chrome, toggle ON, then return to Files and tap the APK again.

"PLAY PROTECT BLOCKED THE INSTALL"

Tap **Install anyway** in the same dialog. If that option is missing, tap the three-dot menu in the corner of the Play Protect screen and choose **Install without scanning**.

"I FORGOT MY PASSWORD RIGHT AFTER CREATING IT"

On the sign-in screen, tap **Forgot Password?** — the app sends a reset link to the email the patient used.

"MY PHONE IS TOO OLD"

If Android version is below 8.0, the install path is not supported. Flag your team lead.

"THE APP KEEPS SHOWING THE PDPA SCREEN EVEN AFTER I TAPPED AGREE"

The patient probably tapped the button without ticking the consent box. Tap the checkbox first, then the button.

"THE APP IS ASKING ME TO UPDATE RIGHT AFTER INSTALL"

Normal. Between the APK build and now, the team may have published an over-the-air update. Tap the banner, the app refreshes itself.

When you are done

DONE CHECKLIST

- ✓ App installed and the patient can sign in on their own.
- ✓ Profile filled in (or honestly blank where they did not know).
- ✓ PDPA consent given.
- ✓ “Install from unknown sources” turned OFF.
- ✓ Patient has seen all four tabs and Privacy Settings.
- ✓ Patient knows they can delete their account and data anytime.
- ✓ Patient knows you do not have access to their password.

Thank the patient for their patience and time. They are now set up to use HaemoCare.

Tear-off quick checklist

Print this page. Use it as a single-page reference at the patient’s side.

- 1 Wi-Fi connected, Android 8+, 200 MB free.
 - 2 **Scan the QR code** with the patient’s camera — opens the GitHub releases page.
 - 3 Tap the **.apk** file under **Assets**; tap **Keep** if Chrome warns.
 - 4 Tap the downloaded APK from notifications or Files.
 - 5 Allow install from this source: **Settings** → **toggle ON** → **back**.
 - 6 Tap **Install**; tap **Install anyway** if Play Protect prompts.
 - 7 Open the app. Toggle EN / TH if needed.
 - 8 **Create Account** — patient types their own password.
 - 9 Sign in. Complete the medical profile. Set days between transfusions.
-

10 PDPA consent — read together, tick the box, tap **I Agree & Continue**.

11 Turn **install from unknown sources** OFF.

12 Walk through Passport, Log, Appointments, History, Privacy Settings.

13 Ask the patient to sign in once on their own.

14 Done.
